

Specification Sheet



Product Name: TelotristatEtiprate

Item No.: MED10055

Purity Specification: $\geq 98\%$

Molecular Formula : C₂₇H₂₃ClF₄N₄O₃

CASNumber: 1040526-12-2 **FormulaWeight:** 562.940

Storage: -20°C

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816