

Specification Sheet



Product Name: (\pm)-Cloprostenol sodium salt

Item No.: MED13135

Purity Specification: 98%

Molecular Formula : C₂₂H₂₈CLNAO₆

CASNumber: 55028-72-3

FormulaWeight: 446.897


Storage: Desiccate at -20 °C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816