

# Specification Sheet



**Product Name:** Presapogenin CP4

*Item No.:* MED16093

**Purity Specification:** 98%

**Molecular Formula :** C<sub>46</sub>H<sub>74</sub>O<sub>15</sub>

**CASNumber:** 75799-18-7

**FormulaWeight:** 867.071


**Storage:** -20 ° C

**Test**

**Specification**

	<b>Min.</b>	<b>Max.</b>
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web:www.med-bio.cn Tel:400-186-1816