

Specification Sheet



Product Name: Arachidic Acid-d3

Item No.: MED26426

Purity Specification: $\geq 98\%$

Molecular Formula : C₂₀H₃₇D₃O₂

CASNumber: 202480-70-4 **FormulaWeight:**

Storage: -20

Test	Specification	
	Min.	Max.
PURITY	98%	100%

Structure :

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816