

Specification Sheet



Product Name: Tangeretin

Item No.: MED80112

Purity Specification: 98%

Molecular Formula : C₂₀H₂₀O₇

CAS Number: 481-53-8

Formula Weight: 372.369


Storage: 2~8 ° C

Test

Specification

	Min.	Max.
PURITY	98%	100%

Structure :



Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Inhibitors / Regulators / Agonists / Intermediates / Plant extracts

Web: www.med-bio.cn Tel: 400-186-1816